BT223

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BT223

Self-Contained Biological Indicator.

Super Rapid Readout Fluorescence System.

Usage

Monitoring gravity air-displacement Steam Sterilization processes at 132 °C – 135 °C.

Applicable Regulation

ISO 11138-1:2006 and ISO 11138-3:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

Characteristics

Polypropylene random copolymer tube: 50.4 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Polypropylene cap: 10.7 mm external diameter, 16.4 mm high. Wall thickness: 0.9 mm

Cap filter: medical grade paper, 60 g/m², 17.0 mm diameter.

Glass ampoule: 35.0 - 40.0 mm high. External diameter: 6.8 mm. Wall thickness: 0.2 – 0.3 mm.

Culture medium 0.5 – 0.7 ml, purple color.

Polypropylene microfiber on spore carrier, 17.0 mm diameter, 80 g/m².

Spores carrier: Filter Paper: 5.0 x 25.0 mm, 160 g/m²

≥ 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Final fluorescence reading is performed after 30 minute-incubation at 60 °C (sensitivity: 97 %).

An optional visual pH color change confirmation could be made after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 30 minute-reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: if 7-day readout is performed, a humidified environment will be required to avoid medium to dry out.

D-Value: Not lower than 1.5 minutes at 121 °C. Another D-value is informed at 132°C.

Z-Value: Not lower than 6 °C.

Environmental conditions during manufacture

Room temperature 15-30°C, relative humidity 30-80%. Sterility conditions are only necessary during the inoculation and manufacturing process, which is performed in laminar flow.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place in its original box.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

24 months.

Packing

10 units per aluminum foil bag. Weight: 35.6 g.

50 units per box. Weight: 142.2 g.

Packing information: product code and description, process for intended use, presentation, classification according to regulation, manufacturer information and data on pack's label.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label. 1.5 mm x 10.0 mm chemical indicator line, printed with Steam reactive ink (Color change to brown). Graph showing final fluorescence reading time, product code, batch number, process for intended use and bacterial strain printed in black.

On product's pack: product code and description, batch number, bacterial strain and load, manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

Incubation at 60 °C must be performed in Bionova® IC10/20FR, IC10/2FRLC or MiniBio Auto-Reader incubator.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.