# BT50/6

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# BT50/6

## Spore Strip Biological Indicator

For Steam and Formaldehyde sterilization processes.

### Usage

Monitoring 121 °C - 134 °C Steam and Formaldehyde Sterilization Processes.

### Applicable Regulation

ISO 11138-1:2006, ISO 11138-3:2006 and ISO 11138-5:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999.

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Characteristics

25 mm x 70 mm non-absorbent material package. Brown code.

5 mm x 25 mm filter paper strip.

$10^6$  *Geobacillus stearothermophilus* ATCC 7953 spores per strip.

Upon completion of the sterilization cycle, the spore strip should be tested in a suitable culture medium at 60 °C. We suggest the use of BIONOVA® MC20 or MC1020 culture media to final reading in 24 hours for Steam sterilization and 48 hours for Formaldehyde sterilization. For conventional culture media like TSB, incubation for 7 days at 60 °C is recommended.

D-Value<sub>STEAM</sub>: not lower than 1.5 minutes. Conditions: 121 °C.

Z-Value<sub>STEAM</sub>: not lower than 6 °C.

D-Value<sub>FORM</sub>: not lower than 6 min. Conditions: 60 °C, 1 mol/l formaldehyde.

### Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

T = 10-30 °C, RH 30-80 %.

### Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### Shelf-life

24 months.

**Packing**

10 units per aluminium foil bag. Weight: 16.1 g.

100 units per aluminium foil bag. Weight: 36.3 g.

**Labelling**

Code and description of the product, batch, expiration date, strain and population of bacteria, presentation, regulation, storage conditions and manufacturer information.

**Possible target markets**

Healthcare and Industry.

**Other important information (optional)**

It is advisable to incubate at 60 °C in BIONOVA® IC10/20 incubator.

Read product's directions for use thoroughly before use.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to radiation, H<sub>2</sub>O<sub>2</sub> or any sterilization process other than Steam and/or Formaldehyde.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.