# BT93/6

## Cupones con Esporas

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# BT93/6

## Spore Strips

For Plasma or Vaporized Hydrogen Peroxide sterilization processes

## Usage

Monitoring Plasma or Vaporized Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) Sterilization Processes.

## Applicable Regulation

ISO 11138-1:2006 and IRAM 37102-1:1999.

## Classification

Class 1, according to risk.

## Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

## Characteristics

25 mm x 70 mm non-absorptive material package. Blue code. Perforated stainless steel coupon 34.0 mm x 7.0 mm x 0.8 mm.

10<sup>6</sup> *Geobacillus stearothermophilus* ATCC 7953 spores per coupon. Upon completion of the sterilization cycle, the spore coupon should be tested in a suitable culture medium at 60 °C. We suggest the use of Bionova® MC20 or MC1020 culture media to final reading in 24 hours. For conventional culture media like TSB, is recommended incubation for 7 days at 60 °C. D-value: at 50 °C, 2mg/l H<sub>2</sub>O<sub>2</sub>.

## Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80%. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

## Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

## Transportation conditions

Storage conditions should be strictly followed. Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

## Shelf-life

24 months.

### **Packing**

100 units per box.

Packing information: product code and description, process for intended use, classification according to regulation, manufacturer information and data on box 's label.

Information on product 's package: product description, process for intended use, batch number, manufacture date and expiration date.

### **Labelling**

On product's box: batch number, manufacture date and expiration date.

### **Possible target markets**

Healthcare and Industry.

### **Other important information (optional)**

It is advisable to incubate at 60 °C in Bionova® IC10/20 incubator.

Read product's directions for use thoroughly before use.

NOTE: When necessary, limits of each product's parameters should be specified.

### **Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to Radiation, EO or any sterilization processes other than H<sub>2</sub>O<sub>2</sub>.

The positive biological indicator can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

Note: The product can be incubated in Bionova® MC20 or MC1020 culture media or other culture medium made up of soybean casein and yeast extract.