# BT95

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# BT95

## Self-Contained Biological Indicator.

*Rapid Readout Fluorescence System.*

### Usage

Monitoring plasma or vaporized Hydrogen Peroxide sterilization processes (H<sub>2</sub>O<sub>2</sub>).

### Applicable Regulation

ISO 11138-1:2006 and IRAM 37102-1:1999.

### Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Characteristics

Polypropylene tube: 50.4 mm high x 8.5 mm external diameter. Wall thickness: 0.5 mm

Polypropylene cap: 10.7 mm external diameter, 16.4 mm high. Wall thickness: 1.0 mm

Cap filter: polyethylene fibers, 17.0 mm diameter.

Glass ampoule: 35.0 – 40.0 mm high. mm. External diameter: 6.8 mm. Wall thickness: 0.2 – 0.3 mm.

Culture medium 0.5 – 0.7 ml, purple color.

Polypropylene microfiber on spore carrier, 17.0 mm diameter, 60 g/m<sup>2</sup>.

Spores carrier: polyethylene fibers: 5.0 x 40.0 mm.

≥ 10<sup>6</sup> *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Final fluorescence reading is performed after 2-hour incubation at 60 °C (sensitivity ≥ 97 %).

An optional visual pH color change confirmation could be performed after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 2-hour reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

D-Value: at 50 °C, 2 mg/l H<sub>2</sub>O<sub>2</sub>.

### Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during inoculation process performed in laminar flow.

### Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place in its original box.

### **Transportation conditions**

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

### **Shelf-life**

24 months.

### **Packing**

10 units per bag. Weight: 35.7 g.

50 units per box. Weight: 149.0 g.

Packing information: product code and description, process for intended use, presentation, regulation, manufacturer information and data on pack's label.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

### **Labelling**

On product: 17.0 mm x 33.0 mm polypropylene label. 1.5 mm chemical indicator line, printed with Hydrogen Peroxide reactive ink (Color change to green). Graph showing final fluorescence reading time, product code, batch number, expiration date, process for intended use and bacterial strain printed in black.

On product's pack: product code and description, batch number, bacterial strain and load, manufacture and expiration date, bar code and datamatrix code.

### **Possible target markets**

Healthcare and Industry.

### **Other important information**

It is advisable to incubate at 60 °C in Bionova® IC10/20FR, IC10/20FRLCD or MiniBio Auto-Reader incubator

Read product's directions for use thoroughly before use.

### **Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry heat, Radiation or any sterilization process other than Hydrogen Peroxide.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer

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