Cleaning Performance Test

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Results reference guide for Ultrasonic Washers

Grade	Visual results
0	Processed and correct
1	
2	
3	
4	
5	Unprocessed

NOTE: These results correspond to examples of possible results when using either CDWA3 or CDWA4 with the corresponding holder for ultrasonic cycles, CDWAH-U. The examples shown are only for approximations for user guidance.

Exclusive use for professionals and Health Care Institutions

Cleaning Indicators

For cleaning performance tests

Control of cleaning and washing processes of medical instruments and materials is highly important since the outcome affects subsequent disinfection and/or sterilization processes success. Cleaning control is particularly relevant in automatic washing machines, where even though process parameters may be well controlled, the cleaning final result depends largely on uncontrolled factors. Washing efficacy is influenced by detergent, water temperature, time and washing penetration. A good cleaning result is achieved by an appropriate combination of these variables. The same considerations apply for washing processes in ultrasonic washing machines.

Product description

Chemdye® CDWA indicators consist of a 70 x 30 mm temperature stable support containing a mixture of specially combined natural colored components (proteins, carbohydrates, fatty acids, and dyes). In all cases, indicator colored test soil is non-toxic, dissolves with most of the cleaning agent and is swept along with the washing fluid. CDWA cleaning indicators formula allows detection of all the factors that affect the outcome of the washing process. CDWA indicators can be used for routine control of cleaning processes and allow monitoring of washing effectiveness in automatic machines including ultrasonic washing machines with different washing configurations. CDWA cleaning system offers two different levels of cleaning with different adhesion characteristics, thus allowing different challenge levels to test cleaning efficacy. The indicators efficiency is validated by comparing against the test soils detailed in ISO/TS 15883-5:2005 standard annexes.

Indicators Classification

Code	Sample color	Annex of standard	Challenge level	Recommended use in thermodisinfection washing	Recommended use in ultrasonic washing (35-40 kHz) (*)
CDWA3	Blue	G	High	Temperature > 40 °C Cycle duration > 3 min.	Temp. 30-50 °C Enz (*) 0.5% ≥ 10 min. Alk (*) 0.2% ≥ 5 min.
CDWA4	Fuchsia	Q	Very high	Temperature > 50 °C Cycle duration > 5 min.	Temp. 30-50 °C Enz (*) 0.5% ≥ 20 min. Alk (*) 0.2% ≥ 10 min.

- *1: The tests carried out with CDWA indicators in ultrasonic washing machines were performed using the holder designed for this purpose, CDWAH-U, which allows to fasten the indicator in a specific location inside the ultrasonic bath, without representing a challenge to the washing processes.
- *2: Enzymatic detergent, tested at a standard concentration of 0.5% V / V.
- *3: Alkaline detergent, tested at a standard concentration of 0.2% V / V.
- *4: Products were not tested with enzymatic detergent in the range between 50-70 °C, because most enzymatic detergents may not work efficiently at this temperature range.

NOTE: It is recommendable to test the 2 types of CDWA indicators along with the instrument, in order to select the one which is more appropriate for the washing machine performance. It is advisable to start the washing control validation process using the most rigorous CDWA4 product. The indicator that washes off completely during the washing/disinfection program should be used for routine monitoring. It should be noticed that the choice of indicator will also depend on the hardness of dirt adhered to the object to be washed. Therefore, the CDWA indicator selected for routine monitoring of a particular washing machine could be changed according to the hardness of the soil on the load. Once the CDWA indicator appropriate for the washing requirements has been selected, we recommend placing several identical indicators at different positions of the washing machine, so as to monitor the performance of all areas of the machine during the washing process.

Characteristics

CDWA indicators have a space for recording information and are placed within a special holder (CDWAH or CDWAH-U holder, see figure) that can be fixed on a tray. Holder design mimics a realistic challenge of thoroughly cleaning hard-to-reach areas to the washing/disinfecting process. When CDWA is placed within the CDWAH holder, half of the test soil is obscured (COVERED ZONE in the figure), thus simulating the joint of a surgical instrument, with water stream hitting it indirectly. The other half of the test soil is covered by a mesh (BASKET ZONE in the figure), thus posing a difficulty for water stream to reach it, mimicking shadowed or occluded surfaces. CDWA indicators provide constant cleaning characteristics compared to the use of natural samples based assays which might suffer modifications of adhesion and polarity properties during storage time. The CDWA system provides high confidence and high reliability to reveal defects in the cleaning process.

Instructions for use

1. Choose the CDWA indicator according to the desired difficulty.
2. Place the indicator on the CDWAH or CDWAH-U holder.
3. Place the system on the tray along with the material.
4. Once the cleaning cycle is over, remove the indicator.
5. Analyze CDWA indicators using the reference table and record results. If unfavorable results are obtained, perform the appropriate corrective actions and repeat the washing process.

Use at least one CDWA indicator for each cleaning cycle.

NOTE: For the holder model in which the test soil is completely covered by a mesh, the effect of washing/disinfecting process is homogeneous, unlike to the latest holder design with mesh only in the top half, which allows appreciating two different effects. Therefore, if the full mesh holder is used, it should be noted that the results shown in the reference table for the top half of the test soil will be repeated in the bottom half.

Observations

Technical specifications outlined in ISO/TS 15883-5:2005 standards describe different types of test samples having different properties and hardness facing cleaning, but no recommendation is made about which sample type should be used. Currently, none of the samples is suggested as a reference or standard because no method has been defined to compare the properties of the materials used in these tests. For this reason, Terragene® provides two washing indicators with different levels of response and adhesion.

Precautions

Do not reuse CDWA indicators. Do not touch the indicator ink with your hands.

Storage

CDWA indicators should be stored at temperatures between 10-30 °C, protected from light and moisture, within its original packaging.
NOTE: in case of remaining traces of indicator colored spot after the washing process, do not store the used indicator under direct light since this may deteriorate its coloring.

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NOTE: This reference table corresponds to results when using CDWA4 in thermodisinfection washers and could be a consequence of failure in one or a combination of variables. When detecting the failed variable/s, make the corrective actions. If the problem persists, contact the technical service of your washing machine. The higher the quantity of X, the higher the influence of that variable on the failure of the cycle. Same results could be obtained when using CDWA3 under the user's specific conditions. Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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